


Human factors lab

Lab location E1025



Device name : Jackson Strength Evaluation System

Used for : to measure the maximum voluntary lifting capacity .

Experiment associated with it :-

The maximum voluntary lifting capacity .



Device name : The Two-Arm Coordination Tester

Used for : The ability to move both arms in a simultaneous and coordinated manner is measured

Experiment associated with it :-

Assesment of dexterity & coordination



Device name :- Minnesota Manual Dexterity Board

Used for : It's used in the manual dexterity test

Experiment associated with it :-

Assesment of dexterity & coordination







Device name : Automatic Tally Maze

Used for : The human learning curve effect can be studied using the tally maze

Experiment associated with it :-

Assesment of dexterity & coordination

	<p>Device name :-The Hole Type Steadiness Tester</p> <p>Used for : It measures the arm hand steadiness</p> <p>Experiment associated with it :- Assesment of dexterity & coordination</p>
	<p>Device name : <u>The Mirror tracer</u></p> <p>Used for : It measures the bilateral transfer of training .</p> <p>Experiment associated with it :- Assesment of dexterity & coordination</p>
	<p>Device name :- Grooved Peg Board</p> <p>Used for : Testing the finger dexterity</p> <p>Experiment associated with it :- Assesment of dexterity & coordination</p>
	<p>Device name : <u>Silent Response Counter</u></p> <p>Used for : It's used to measure the number of errors in dexterity and coordination test exp.</p> <p>Experiment associated with it :- Assesment of dexterity & coordination</p>

	<p>Device name :- Hand Grip Dynamometer</p> <p>Used for : Measuring the maximum voluntary contraction of the biceps muscles in the arm .</p> <p>Experiment associated with it :- maximum voluntary contraction (MVC)</p>
	<p>Device name :_Single Event Stop Watch</p> <p>Used for ;- To measure the time</p> <p>Experiment associated with it :- Assesment of dexterity & coordination</p>
	<p>Device name :- Professional Vision Tester</p> <p>Used for : It's used to screen different vision defects</p> <p>Experiment associated with it :- Vision tester</p>
	<p>Device name :_ Multi-choice reaction time apparatus</p> <p>Used for: to measure the response time to different stimuli.</p> <p>Experiment associated with it :- Response time</p>



Device name :- Anthropometry Set

Used for : It's used to take the required body dimensions which are used in the redesign process

Experiment associated with it :-
Anthropometry and work space design



Device name :_Television & video

Used for : It's used to watch video tapes about special topics in the ergonomics field